Change Log for Abstract Plus Installations

Version 1.14.18 and Prior

Layout: NAACCR 10.2 and prior

Release Date: October 21, 2005 and prior

Changes from Previous: Layout

Version 2.0.11

Layout: NAACCR 10.2 Release Date: Beta Changes from Previous:

- Automatic Internet upgrade option
- Improved menus with additional options
- Microsoft[®] Access 2000 compatible
- Configuration selection sets display type, edit set, and Collaborative Staging (CS) calculation
- Configuration remains active until another is selected
- Added the ability in Preferences to allow directory selection for database, backup, and restore
- · Able to create sections, which are used to group display fields
- · Expanded exporting options
- Expanded print options
- Expanded keyboard navigation
- Added site-specific CS and surgery code lookups
- Text fields available for critical field selection in Project definition
- Added the ability to sort field title and pop-up lookups
- Added new descriptive messages for fields in abstract grid without a choice list
- Doctor Query (DQS) available without open abstract or physician field selection
- Added the ability to set the default value of null for a protected field
- Added Vital Status to selected cases report and abstract selection
- Added the options in exporting a file, to export by DxDate, view and print selections before export, and continue or reselect records
- Added the ability in Preferences to allow the choice of NAACCR or delimited, which allows selection of fields and format when exporting
- Added the ability in Preferences to allow F5 key selection to copy Dx address fields to Curr address fields
- Added the ability in Preferences to allow personal identifier (PID) instead of social security number (SSN) and language selection for international use
- Added the ability in Preferences to allow only the fields in the display type to show in the EDITS

Version 2.1.4

Layout: NAACCR 11.0

Release Date: February 16, 2006

Changes from Previous: Enhanced features with NAACCR11.0 layout.

Version 2.1.5

Layout: NAACCR 11.0

Release Date: February 28, 2006

Changes from Previous: Exe changed to handle blank NAACCRRecVersion; not included in display

with default value.

Version 2.1.6

Layout: NAACCR 11.0

Release Date: March 10, 2006

Changes from Previous: Exe changed to include a message when F5 pressed for derived collaborative

staging with a blank or incorrect DxDate.

Version 2.1.7

Layout: NAACCR 11.0

Release Date: March 16, 2006

Changes from Previous: Exe changed for password-protected database issues.

Version 2.1.8

Layout: NAACCR 11.0

Release Date: March 28, 2006

Changes from Previous: Exe correction for configuration change with different edit set; correction for

completed status retained instead of changing all to incomplete.

Version 2.1.8

Layout: NAACCR 11.0 Release Date: April 17, 2006

Changes from Previous: Include NAACCR11.0B EDITS metafile.

Version 2.1.9

Layout: NAACCR 11.0 Release Date: May 11, 2006 Changes from Previous:

Exe changed to include:

- Message correction for use of removable drive for backup
- ProjID accepts 01–09, changes to correct form of 1–9
- Current edit set included in message box text for abstract with errors
- Added 'Held by Abstractor' field to open/find abstract list
- Added sort preference capability for pick-list tables; default is sorted by description text if list greater than 13 items, except for races, otherwise sorted by code.

Master.mdb changed to correct error with 'All Fields' display.

Version 2.1.10

Layout: NAACCR 11.0 Release Date: May 16, 2006

Changes from Previous: Exe changed to not correct for invalid day, but pass to EDITS for error

message.

Version 2.1.11

Layout: NAACCR 11.0 Release Date: July 21, 2006 Changes from Previous:

Page 2 of 8 Updated: November 19, 2009

Exe updated to include:

- Correction to date editing: certain months with unknown day and known year would return an error, along with the year 1999, century now used in finding unknown year.
- Correction for on export, project ID recorded in Hormone Text field. On import, project ID was assigned as 1 by default or selected values, so value never lost.
- Protected database upgrade correction.
- Upgrading message changed for correct version.

Version 2.1.12

Layout: NAACCR 11.0

Release Date: August 23, 2006 Changes from Previous:

Exe updated to include:

- For database folder selection, message for drive not available, allowing a new selection.
- Alphanumeric allowed for hospcode when edits allows.
- AutoFillAbstractor option applied only on new abstract.
- When state = XX, allow entry of 3-digit country code in county field.
- RxDateSystemic, RxDateRad; any date field allowed 88888888, saved.
- Invalid histology/psite gives error message for CStage lookup.

Version 2.1.13

Layout: NAACCR 11.0

Release Date: August 31, 2006 Changes from Previous:

Exe updated to include:

- Correction for blank vital status to list as blank in reports.
- Correction for alphanumeric hospcode guerying.

Version 2.1.14

Layout: NAACCR 11.0C

Release Date: September 27, 2006

Changes from Previous:

Include NAACCR11.0C EDITS metafile, ver 010300 cstage.dll, and csdllvb.dll

Exe updated to include:

- Correction for EDITS return on protected field, in abstract.
- Correction for selection of Cancel button on backup prompt.
- Error message not generated for F5 calculation of DxAge.
- View of large export file allowed in export log.

Rpsetup_Manual.exe – external program that can be copied to the application folder with an upgrade.zip file containing replacement files, that unzips and copies files in upgrade.zip to appropriate folders. A manual upgrade that does not require administrative privileges as uninstall/install, but requires current version of at least 2.1.4.

Version 2.1.15

Layout: NAACCR 11.0C

Release Date: October 27, 2006

Changes from Previous:

Exe updated to include:

 Creates a separate table to save preferences for database folder, to preserve selections with an upgrade.

Version 2.1.16

Layout: NAACCR 11.0C

Release Date: November 3, 2006

Changes from Previous: Exe correction needed for layout upgrade to NAACCR11 only.

Version 2.1.17

Layout: NAACCR 11.0C

Release Date: November 13, 2006

Changes from Previous:

Exe updated to include:

Correction for backup to external drive.

- Enhancement to export record type recorded in abstract, if other than 'A'.
- Correction to not allow deletion of only userID.
- Correction to allow "mm/dd/yyyy" format for default date setting.

Version 2.1.18

Layout: NAACCR 11.0C

Release Date: December 11, 2006

Changes from Previous:

Exe updated to include:

- Message given on FUPContNumandSt field for F5 copy of Dx address fields, like DxCurrNumandSt, if allowed by options.
- On import of abstract, if DateExported exist then DateOffLoaded marked.
- Allow import of special "M" record type, as normal "A" record type.
- Correction for export log records without SQLtext field (prior to v2.1);message given instead of runtime error.
- Correction for saving invalid values for derived CStage fields calculated when not displayed.

Generic install database changed to include a configuration with all CStage calculations.

Version 2.1.19

Layout: NAACCR 11.0C

Release Date: December 15, 2006

Changes from Previous:

Exe updated to include:

• Correction for importing, exporting, or saving text fields with special characters from copy/paste from an outside source that caused multi-line records for an abstract on the export text file.

Version 2.1.20

Layout: NAACCR 11.0C

Release Date: January 10, 2007

Page 4 of 8 Updated: November 19, 2009

Changes from Previous:

Exe updated to include:

- Check box option on import for 'on import of abstract, if DateExported exist then DateOffLoaded marked'. Not checked by default.
- Correction to clear previous export log record details, when search by abstract.

Version 2.1.20

Layout: NAACCR 11.0C

Release Date: January 18, 2007

Changes from Previous:

Exe updated to include:

 Request for incomplete status on import as default value (previously complete by default), selection available; RPUG did not disagree, if needed preference setting can be created.

Version 2.1.21

Layout: NAACCR 11.0C Release Date: March 6, 2007 Changes from Previous:

Exe updated to include:

- Correction for site-specific surgery pick list headings and notes omitted; codes not changed.
- CStage sub-titles added for site-specific factors.
- Correction to allow edit of fields added up to the end of the record.

Master update:

• Typo corrected in Surgery_Bones table: 400.0 to 40.0.

Version 2.1.22

Layout: NAACCR 11.0C Release Date: March 14, 2007 Changes from Previous:

Exe updated to include:

• Correction for EDITS not run initially on Abstracted by field for new abstract.

Version 2.2.1

Layout: NAACCR 11.1 Release Date: April 19, 2007 Changes from Previous:

Exe updated to include:

- Export report can be sorted by column selection in print preview.
- Option to query a group of abstracts for deletion, and deselect individuals from the group. Backup of database created that can be restored until next deletion.

Updated online help for NAACCR v11.1.

Updated nconv.dll (v4.2.0) for NAACCR v11.1.

Version 2.2.1

Layout: NAACCR 11.1 Release Date: May 10, 2007

Changes from Previous: Updated Edits engine (v2.0.2.0).

Version 2.2.2

Layout: NAACCR 11.1 Release Date: June 20, 2007 Changes from Previous:

Exe changes:

- Correction for Count of Cases Summary by month when 'All cases' selected, with several years data. Count was correct when date range used.
- Correction for EDITS, Re-Edit, and "Batch Completion Re-check..." error count, made the same for all instances. SuppressInit option was not always applied. SuppressInit default changed to False, which allows all edit errors to display, even if all fields in edit error not in display.
- Allow removal of invisible or protected property for display fields, through administration screens.

Version 2.2.3

Layout: NAACCR 11.1 Release Date: July 24, 2007 Changes from Previous:

Updated edits metafile (v11.1B) – correction for wrong edit description displayed with use of v11.1A and EditWriter3.

Exe changes:

- Correction for runtime error for F4 surgery lookup for Pancreas Psite only.
- Message box added to confirm doctor name for DQS selection. If selection click not in DocID column then the first DocID is selected.
- Changes to Delete Abstract screen, including DateOffLoaded added to selection criteria, toggle of the select column heading to change each to opposite selection, and final message box confirming abs IDs to delete.
- "Note:" in select value box changed to clarify that this is a message following and not a pick-list (Note: No Choice List for this field, but some available values to enter or additional helpful hints.). Also changed to bold font and yellow background for distinction.
- Correction for text fields saving with page up/down and home/end keys.

Version 2.2.4

Layout: NAACCR 11.1

Release Date: October 2, 2007 Changes from Previous:

Exe changes:

- Correction for cutoff of text of some CStage pick list rows (ie. C569 CSMetsDX; C221 CSLymphNodes; C329- CSSSF5)
- Option for titlebar name added to DBOptions table; allows icon/titlebar to show customized Abstract Plus name for multiple installs on workstation; Installation options allow for multiple installs with separate Start menus.

Page 6 of 8 Updated: November 19, 2009

Correction for incorrect query selection of abstracts on Export Report if re-export of abstracts
previously marked for export and the check for mark for abstract selected. Correct files were
exported with incorrect report listing. Also corrected concatenation of export type for log of export.

Version 2.3.0

Layout: NAACCR 11.2

Release Date: February 1, 2008

Changes from Previous:

Includes NAACCR11.2 EDITS metafile and ver 010400 cstage.dll.

Exe changes:

- Correction for pathname in backup prompt with selected database folder.
- Export Abstracts, Selected Cases Report, and Print Range of Abstracts options changed to allow entry of a range of Abstract IDs for non-consecutive numbers, up to 20 characters (ie. 2,4,6–10, 13–15, 20).

Version 2.3.1

Layout: NAACCR 11.2 Release Date: July 14, 2008 Changes from Previous:

Includes version 010401 cstage.dll and version 11.2 online help.

Exe change:

• Correction for export by projectID.

Version 2.4.0

Layout: NAACCR 11.3 Release Date: March 5, 2009 Changes from Previous:

Abstract Plus upgrade for NAACCR11.3

Version 2.4.1

Layout: NAACCR 11.3

Release Date: November 9, 2009

Changes from Previous:

Includes updated Microsoft files.

Exe change:

Correction for display type ability to add any field.

Version 2.4.2

Layout: NAACCR 11.3

Release Date: November 19, 2009

Changes from Previous:

Exe change:

Correction for display of wrong surgery primary site lookup table reported for C171/9930.