



# MMWR™

Morbidity and Mortality Weekly Report

**MMWR Dispatch**  
**Vol. 53 / March 23, 2004**

## Notice to Readers

### **Update: Manufacturer's Recall of Rapid Cartridge Assay Kits on the Basis of False-Positive *Cryptosporidium* Antigen Tests**

On March 4, 2004, CDC announced that a manufacturer had voluntarily recalled rapid cartridge assay kits because of false-positive *Cryptosporidium* antigen tests (1). An additional lot of a *Cryptosporidium*/*Giardia* rapid assay has been recalled voluntarily from laboratories by the distributor (Meridian Bioscience, Inc., Cincinnati, Ohio) on the basis of their findings that *Cryptosporidium*-negative samples were weakly reactive with this lot (ImmunoCard STAT!®, lot no. 081138 [expires October 5, 2004]). CDC recommends reconfirmation of positive test results (by using direct fluorescent antibody testing or modified acid-fast stained smears) obtained with ImmunoCard STAT!® rapid assays from all recalled lots.

#### **Reference**

1. CDC. Manufacturer's recall of rapid cartridge assay kits on the basis of false-positive *Cryptosporidium* antigen tests—Colorado, 2004. *MMWR* 2004;53:198.

All *MMWR* references are available on the Internet at <http://www.cdc.gov/mmwr>. Use the search function to find specific articles.

Use of trade names and commercial sources is for identification only and does not imply endorsement by the U.S. Department of Health and Human Services.

References to non-CDC sites on the Internet are provided as a service to *MMWR* readers and do not constitute or imply endorsement of these organizations or their programs by CDC or the U.S. Department of Health and Human Services. CDC is not responsible for the content of these sites. URL addresses listed in *MMWR* were current as of the date of publication.