

## Update Instructions for MMDS 2007.02 ACME/TRANSAX Update

This program update is being distributed using the internet, instead of on CD. To update this program, your PC will need to be internet capable or you will need to copy the files to a disk.

- **If you do not do ACME coding, you do not need to run this update.**
- **Before beginning this update make sure ACME/TRANSAX is not running on your system. If ACME is running, exit that program before performing the update.**
- **Also, this is an update only, not a full install. You will need to have completed the 2007.01 installation before running this update.**

To update the programs, direct your browser to the NCHS website at:

<http://www.cdc.gov/nchs/about/major/dvs/medsof.htm>


You can also open SuperMICAR or any of the MMDS programs and find this link under the HELP menu by clicking on **About**.

Once there, you can click on the UPDATE 2007.02 link and be taken to the 2007.02 folder where you will see the file UPD-0703.EXE, as shown below.

Double-click on the file to start the update..

Depending on your system you may get the file download box shown below. Just click on Open if you do.

The PKZIP box will appear as shown below.



The default location is shown under the “Extract to” label and is C:\MMDS2007. You can change this by using the menus to the right and below that window. When you have the location you want, click on the Extract button in the lower left to complete the update.

It should take no more than a second or two for the program to update. The PKZIP box will vanish when completed.

To be sure everything went well, click on your ACME 2007 desktop icon and verify that the version is now 2007.02.